

WECK



Weck

Hem-o-lok Ligation System

Safety Without Compromise



Teleflex



The Clip Bank

Easy Loading with Convenient Adhesive Strip

- Our Hem-o-lok ligation system combines speed and safety in one single procedure.
- Choose from a variety of ligating clip sizes for the ligation of vessels and/or tissues between 2 and 16 mm.
- The Hem-o-lok clip bank allows for effortless and precise loading, thus making surgical interventions more efficient.
- Closure of the clip is confirmed by tactile and audible feedback.
- The colour code of the clip bank corresponds to that of the applier.
- Each cartridge contains 6 clips.
- A special adhesive strip at the bottom of each cartridge enables you to choose the best possible site for clip delivery.
- Unlike absorbable clips, Hem-o-lok ligating clips can be removed with our specially developed removers if necessary.
- The clips are inert, non-toxic and radiolucent.
- The cartridges adhere to any sterile surface.

The Appliers

Performance Features and Cleaning Options

- Thanks to the special clip design, closure of the clip can be clearly felt by the surgeon during application.
- Both curved and angled appliers are available for open surgery.
- Endoscopic Hem-o-lok appliers are available in sizes ML, L, and XL. They can be taken apart for easy cleaning, thus eliminating potential contamination.
- A 360° rotation knob facilitates access and improves visibility.
- The appliers are colour-coded to match the colours of the corresponding cartridges.

Ligating Clips



Open Surgery



Hem-o-lok Ligating Clips

REF.	CLIP SIZE	LIGATION OF VESSELS AND TISSUE	CLIPS/CARTRIDGE	CARTRIDGES/BOX	CLIPS/BOX
WK544230	● ML	3–10 mm	6	14	84
WK544240	● L	5–13 mm	6	14	84
WK544250	● XL	7–16 mm	6	14	84

Appliers

REF.	CLIP SIZE	WORKING LENGTH	QTY
Curved			
WK544170	● ML	20 cm	1
WK544180	● L	20 cm	1

Appliers, Long

REF.	CLIP SIZE	WORKING LENGTH	QTY
Curved			
WK544171	● ML	27 cm	1
WK544181	● L	27 cm	1
WK544191	● XL	27 cm	1

Appliers, Long, 70° Angled

REF.	CLIP SIZE	WORKING LENGTH	QTY
WK544172	● ML	27 cm	1
WK544179	● L	27 cm	1
WK544192	● XL	27 cm	1

Clip Removers

REF.	CLIP SIZE	WORKING LENGTH	QTY
WK544123	● M ● ML ● L	27 cm	1
WK544124	● L ● XL	27 cm	1

● M medium ● ML medium-large ● L large ● XL extra-large

All Hem-o-lok ligating clips are sterile packaged. Ligating clip systems differ in their clip-closure characteristics. It is the responsibility of the user to confirm security of clips on vessels or other tissue. Not intended for contraceptive tubal occlusion.





Hem-o-lok Non-Take-Apart Applier, Long

REF.	CLIP SIZE	Ø APPLIER	Ø TROCAR	WORKING LENGTH	QTY
Straight					
WK544965L	● ML	5 mm	5 mm	45 cm	1
WK544995	● L	10 mm	10 mm	45 cm	1
WK544990	● XL	10 mm	10 or 12 mm	45 cm	1
Curved					
WK544995L	● L	10 mm	10 mm	45 cm	1
WK544990L	● XL	10 mm	10 or 12 mm	45 cm	1



Hem-o-lok Auto Endo 5 Appliers

REF.	CLIP SIZE	Ø APPLIER	Ø TROCAR	WORKING LENGTH	QTY
Straight					
WK543965	● ML	5 mm	5 mm	33 cm	3



Clip Removers

REF.	CLIP SIZE	Ø REMOVER	Ø TROCAR	Ø SILS	WORKING LENGTH	QTY
Straight						
WK544121	● ML	5 mm	5 mm	—	45 cm	1
Curved						
WK544130	● L/XL	10 mm	10 or 12 mm	—	45 cm	1

The Hem-o-lok Ligation System – Precision in Every Detail

Teleflex is a leading manufacturer of metal and polymer clips. Our Weck Hem-o-lok ligation system forms a line of instruments perfectly adapted to each other, meeting all the special needs of surgeons, nurses and operating theatre personnel alike.



The Clip – Top-Class Functionality

- Precision bosses engineered for maximum security.
- Penetrating lock: 1 x 2 teeth to provide effective closure.
- Bowed design enables added tissue to be encompassed.
- Tactile and audible response confirms that the Hem-o-lok ligating clip has locked onto vessel.
- Integrated teeth help prevent slippage in any direction.
- The surface properties of the non-absorbable polymer prevent infection.
- Living hinge provides flexibility when opening and closing.

CLIP SIZE	LIGATION OF VESSELS AND TISSUES
● Medium	2–7 mm
● Medium-large	3–10 mm
● Large	5–13 mm
● Extra-large	7–16 mm

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow, Deknatel, Hudson RCI, LMA, Pilling, Rüsch and Weck – trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Teleflex Medical Private Limited

"Blue Haven" 19 · Harrington Road · Chetpet

Chennai 600031 · India

Tel. +91 (0)44 2225 0340, +91 (0)44 2225 0341, +91 (0)44 2225 0342

www.teleflex.com

Teleflex, the Teleflex logo, Arrow, Deknatel, Hem-o-lok, Hudson RCI, LMA, Pilling, Rüsch and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries.

Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (05/2016). Subject to technical changes without further notice.

© 2016 Teleflex Incorporated. All rights reserved.

94 05 70 - AP 00 ID · REV A · MC/PDF · 05 16 PDF